

medaSEPT NITRILE PREFER BLACK PF

MEDICAL AND PROTECTIVE NITRILE
POWDER-FREE GLOVES IN BLACK



SPECIFICATION	
RAW MATERIAL	NITRYL
COLOR	BLACK
POWDERING AGENT	NONE, POWDER-FREE GLOVES
SHAPE	UNIVERSAL
CUFF	ROLLED
INTERIOR SURFACE	CHLORINATED
EXTERNAL SURFACE	FINGERTIP TEXTURED
AQL	1.0 final product control level: G-1 according to ISO 2859-1

CLASSIFICATION & NORMATIVE STANDARDS		
CLASSIFICATION	MEDICAL DEVICE	PERSONAL PROTECTIVE EQUIPMENT
	CLASS I	CATEGORY III, TYPE B
	in accordance with Regulation (EU) 2017/745	in accordance with Regulation (EU) 2016/425
STANDARDS	EN 455 1-4	EN ISO 374-1, EN ISO 374-2,
	EN ISO 15223-1, EN 1041	EN 16523, EN ISO 374-4,
		EN ISO 374-5, EN ISO 21420
CODE BASIC UDI-DI		590173887RDNPLG

PRODUCT PROPERTIES	
RESISTANCE TO BACTERIA AND FUNGI	in accordance with EN ISO 374-2
RESISTANCE TO VIRUSES	in accordance with EN ISO 374-5 (ISO 16604) & ASTM F 1671
CHEMICAL RESISTANCE	in accordance with EN 16523-1 & EN ISO 374-1
RESISTANCE TO CYTOSTATICS	in accordance with ASTM D 6978

ISO 374-1/TypeB

KPT

ISO 374-5:2016

VIRUS

TECHNICAL PARAMETERS & MASTER DATA						
size		XS (5-6)	S (6-7)	M (7-8)	L (8-9)	XL (9-10)
reference number		MSNPREBPF100XS	MSNPREBPF100S	MSNPREBPF100M	MSNPREBPF100L	MSNPREBPF100XL
hand width	[mm]	≤80	80±10	95±10	110±10	≥110
length	[mm]	min. 240				
thickness (single glove wall)	[mm]	finger: 0.09±0.03 palm: 0.08±0.03 cuff: 0.07±0.03				
breaking force	[N]	before aging: 6.0 N after aging 6.0 N				
powder content	[mg/gloves]	< 2 mg/gloves, powder-free gloves				
latex protein content	[µg/g]	none, latex-free gloves				
quantity of pieces in a box		100 pieces according to weights				
quantity in a carton		10 boxes				
EAN box		5901738874855	5901738874862	5901738874879	5901738874886	5901738874893
EAN carton		5901738874909	5901738874916	5901738874923	5901738874930	5901738874947
box size	[mm]	225 (L) X 115 (W) X 60 (H)				
carton size	[mm]	310 (L) X 240 (W) X 235 (H)				
expiration date		5 years from the production date				

medaSEPT NITRILE PREFER BLACK

medaSEPT[®]
QUALITY TO PROTECT

INTENDED USE

A single-use device intended to be used to protect the patient and user in order to minimize the risk of cross-infection. Gloves suitable for medical examinations and other diagnostic procedures.

STORAGE INSTRUCTIONS

Protect from sunlight, fluorescent light and ozone and store in a dry place. The recommended relative humidity of the room where the gloves are stored is 60±20%. Keep gloves away from heating devices and sources of fire.

UTILIZATION

Dispose of in accordance with applicable regulations.

CONTACT WITH FOOD

Suitable for contact with food in accordance with: Regulation (EC) 10/2011, Regulation (EC) 1935/2004.



PRODUCER

medaSEPT S.A.
street Forteczna 19
61-362 Poznan
www.medasept.com

